

by anandic

EasyOne FlowTube

Breathing Tube for
Infection Control

n d d
new diagnostic design

**Avoid cross-contamination
and reduce cleaning to an
absolute minimum with
the single-patient-use
EasyOne FlowTube.**



All parts exposed to patient's breath
are **single-use** and **disposable**

No expiration date

No special storage conditions

Technical data

Material	EasyOne FlowTube: Polypropylen (PP) Bag: Polyethylen (PE)
Transport and storage information	Temp. -20 to 50°C Rel. Humidity 5 to 90% Athmosph. Pressure 500 to 1060 hPa no expiration date
MDD 93/42/EEC	Class I
Dimension (single EasyOne FlowTube)	21 x 38 x 150 mm (H x W x D), 8 g
Compatible with	EasyOne Air

Order information

Part Number	Description
5050-50	EasyOne FlowTube, individually wrapped, dispenser box of 50 pcs
5050-200	EasyOne FlowTube, individually wrapped, box of 200 pcs (dispenser box included)
5050-500	EasyOne FlowTube, individually wrapped, box of 500 pcs (without dispenser box)