

by anandic

OneTouch **NMT**[™]

STIMPOD
NMS 450X



Simple. Economical. Proven.

20 Years of Quantitative **NMT** Monitoring Innovation

The **STIMPOD** Quantitative **NMT** Monitor

Train of Four | Post Tetanic Count | Double Burst | Tetanus | Twitch | Nerve Mapping & Locating

Simple

- OneTouch NMT™ Monitoring
- Quick Setup & No Calibration
- Electrode Placement Verification
- Handheld or IV Pole Mount
- Data Connectivity
- Versatile Sensor Placement:
thumb, toe or face

Economical

- Reusable 3D AMG Sensor
- Standard ECG Electrode Compatible
- Quantitative NMT Monitoring reduces incidence of residual paralysis, increases operating room throughput and lowers cost of patient care

Proven

- Over 20 years of AMG research in the public domain
- Residual paralysis is identified in 97% of patients
- Trusted by anaesthesiologists in over 40 countries
- STIMPOD technology pioneering quantitative monitoring since 1998



Powered by Proven **AMG Technology**

STIMPOD NMS 450X

Simple. Economical. Proven.

www.xavant.com
sales@xavant.com
+27 12 743 5959

X·AVANT
TECHNOLOGY

FDA
CLEARED
CE
CERTIFIED