

vyaire™  
MEDICAL

by anandic<sup>+</sup>



GreenLight™  
Laryngoscope handles and blades

## Designed to make reprocessing easy

- Stainless steel and plastic ISO standard connection single-patient-use blades reduce time and effort associated with reprocessing of reusable blades
- Standard AA batteries\* (*alkaline or lithium*) used in traditional handle for easy replacement (*single CR-123 lithium battery used in stubby handle*)
- With our fiber optic system there is no longer a need for bulb replacement

## Designed for performance

- Ergonomic and over-molded handles are one of the lightest in the market
- Handle has a durable, slip-resistant grip and a wide, sturdy base
- Bright LED light and no back-flash glare
- Easy to remove battery and light cartridge for quick sterilization between cases. Handle is autoclavable

## Designed to last

- Efficient technology extends battery life beyond that of a standard incandescent bulb light source
- Five-year warranty on all handles



Standard handle with single-patient-use stainless steel blade



Stubby handle with single-patient-use stainless steel blade

Part no.	Description	Handle type	Case qty
4558GSP	GreenLight™ ISO standard laryngoscope handle (Santoprene over-molded handle) - reusable	Standard	1
4559GSP	GreenLight™ ISO standard laryngoscope handle (Santoprene over-molded handle) - reusable	Stubby	1
4559LED	Replacement battery/LED cartridge for GreenLight™ (ISO green standard laryngoscope handle)	Standard	1
4559LEDS	Replacement battery/LED cartridge for GreenLight™ (ISO green standard laryngoscope handle)	Stubby	1



Standard handle



Stubby handle



Stainless steel blades



Plastic blades

Part no.	Description	Blade type	Blade size	Case qty
4602	GreenLight™ ISO standard, stainless steel, SPU laryngoscope blade	Mac	2	20
4603EU	GreenLight™ ISO standard, stainless steel, SPU laryngoscope blade	Mac	3	20
4604EU	GreenLight™ ISO standard, stainless steel, SPU laryngoscope blade	Mac	4	20
4610	GreenLight™ ISO standard, stainless steel, SPU laryngoscope blade	Miller	0	20
4611EU	GreenLight™ ISO standard, stainless steel, SPU laryngoscope blade	Miller	1	20
4612	GreenLight™ ISO standard, stainless steel, SPU laryngoscope blade	Miller	2	20
4613	GreenLight™ ISO standard, stainless steel, SPU laryngoscope blade	Miller	3	20
4002GSP	GreenLight™ ISO standard, plastic, SPU laryngoscope blade	Mac	2	20
4003GSP	GreenLight™ ISO standard, plastic, SPU laryngoscope blade	Mac	3	20
4004GSP	GreenLight™ ISO standard, plastic, SPU laryngoscope blade	Mac	4	20
4010GSP	GreenLight™ ISO standard, plastic, SPU laryngoscope blade	Miller	0	20
4011GSP	GreenLight™ ISO standard, plastic, SPU laryngoscope blade	Miller	1	20
4012GSP	GreenLight™ ISO standard, plastic, SPU laryngoscope blade	Miller	2	20
4013GSP	GreenLight™ ISO standard, plastic, SPU laryngoscope blade	Miller	3	20



For more information, please contact your local Vyairé representative.



**ANANDIC MEDICAL SYSTEMS AG**

Stadtweg 24, 8245 Feuerthalen

[www.anandic.com](http://www.anandic.com)

[info@anandic.com](mailto:info@anandic.com)

Tel. 0848 800 900

**GLOBAL HEADQUARTERS**

Vyairé Medical, Inc.  
26125 N. Riverwoods Blvd.  
Mettawa, IL 60045  
USA



EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands



Vyairé Medical, Inc.  
26125 North Riverwoods Blvd  
Mettawa, IL 60045  
USA



[vyaire.com](http://vyaire.com)

For Global Distribution.

© 2019 Vyairé. Vyairé and the Vyairé logo and all other trademarks or registered trademarks are property of Vyairé Medical, Inc., or one of its affiliates. Medical devices class IIa according to Medical Devices Directive 93/42/EEC. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labelling. VYR-GBL-1900115